


# Vibrasat® Pro: Discover the Power of Gentle Vibrations



Vibrasat® Pro is our latest and most gentle power-assisted liposuction device. The easy-to-use design significantly reduces fatigue, making it ideal for long procedures and helping professionals improve their patients outcomes.



### *Gentle procedure*

Reduced trauma to surrounding tissues

### *Low Noise*

Provides a quiet and stress-free environment for the surgeon

### *QuickLock® System*

Efficient and straightforward cannula change for optimal treatment continuity

### *Vibration/Boost*

3,000 - 5,000 strokes/min at approx. 3 mm  
Boost mode with up to 6,000 strokes/min

### *Anatomically shaped*

Ensures precise fat removal and enhanced surgeon comfort during extended procedures

### *Low heat generation*

during long procedures

### *Ergonomically engineered*

Supporting the surgeon's arm to prevent numbing

### *Fully machine washable and autoclavable*



Follow us on



Amazing what's possible.

# Revolutionize Lipedema Care with Power-Assisted Liposuction



Liposuction can reduce the volume and weight of the affected limbs, improve the shape and appearance of the body, and alleviate pain and discomfort. Power-assisted liposuction is a promising technique that can enhance the safety and efficacy of liposuction for lipedema patients.

Power-assisted liposuction uses high-frequency vibrations to break up and remove stubborn fat deposits and is less invasive than traditional liposuction. When applied to lipedema treatment, the method offers several advantages:

## *Precision*

The vibration of the cannula empowers the surgeon to operate with heightened precision, ensuring the more accurate removal of fat tissue - a critical advantage in areas with delicate anatomical structures.

## *Reduced Physical Exertion*

The innovative technology of Vibrasat® Pro revolutionizes liposuction by significantly reducing the physical strain on surgeons, making the procedure more efficient and comfortable.

## *Shorter Operation Time*

Due to the efficient removal of fat, the operation time can be shortened, which is beneficial for both the surgeon and the patient.

## *Improved Patient Outcomes*

A faster healing process and fewer post-operative complaints lead to more satisfied patients, which in turn can enhance the reputation of the surgeon.

## *Versatility*

Power-assisted liposuction can be used in various areas of the body, making it a versatile option for surgeons offering a wide range of cosmetic procedures.

## *Lower Risk of Irregularities*

Vibrasat® Pro may help reduce the risk of skin irregularities or dimpling post-surgery because the fat is removed more evenly.

